


[illegible]

GPCR-1 (top): Shows the activation of Gq. The ligand binds to the GPCR, causing a conformational change that activates the Gq protein. The activated Gq protein then activates PLC, which in turn activates PKC and Ca²⁺ release.


GPCR-2 (middle): Shows the activation of Gi. The ligand binds to the GPCR, causing a conformational change that activates the Gi protein. The activated Gi protein then activates PI3K, which in turn activates Akt and Rho GTPase.

GPCR-3 (bottom): Shows the activation of G12/13. The ligand binds to the GPCR, causing a conformational change that activates the G12/13 protein. The activated G12/13 protein then activates Rho GTPase, which in turn activates Rac and JNK.

The diagram illustrates the human genome, showing 22 pairs of autosomes and the sex chromosomes (X and Y). Each pair is represented by two homologous chromosomes, one in black and one in color. The chromosomes are arranged in a karyotype format, with pairs numbered 1 through 22, followed by the sex chromosomes (X and Y). The diagram includes labels for each chromosome pair and a legend for the color coding.





Hi-C





Lieberman-Aiden 2010

Hi-C




Lieberman-Aiden 2010

Dovetail Sequencing (Putnam 2015)




10X Genomics








- 60 kb deletion



TTGAGAT
TATGAAGC
TAAATCTC
TACCACCT
GCTGAAB
ATTCCCT
TCTGGGA
GAAATT
TGTTGA
AAGGAG
TTTGGG
CGCCAG
TCCCGC
AATTGC
TCTCC
AAGGCT
AATTGA
GCACAA
ATACCA
GCTTTT
TTTAT

SUMMARY PHASING STRUCTURAL VARIANTS chr2+29989813-chr2+30039813,chr2+42493841-chr2+42543841





MeDIP-Seq DIP-seq

Methylated DNA immunoprecipitation (MeDIP-seq), DNA immunoprecipitation followed by high throughput sequencing (DIP-seq)



hMeDIP-Seq

Hydroxymethylated DNA immunoprecipitation combined with next generation DNA sequencing (hMeDIP-seq)



MDBCap-seq Methyl-Cap-seq MDB-Seq MBDCap-seq MIGS

Methyl-CpG binding domain-based capture and sequencing (MBDCap-seq). Capture of methylated DNA using the MBDCap domain of MeCP1 (MethylCap-seq), MBDCap-seq, MBDCap-seq, MBDCap-seq



BisChIP-Seq ChIP-BS-seq

Bisulfite-treated chromatin immunoprecipitation (BisChIP-seq) and ChIP-seq, to correlate protein modifications with DNA methylation



DNA-Protein Interactions

DNase-Seq DNaseI-Seq

DNase I hypersensitive sites sequencing (DNase-seq)



MAINE-Seq MNase-Seq Nucleo-Seq

Micrococcal nuclease digestion (MNase-seq), Micrococcal nuclease digestion (MNase-seq), Micrococcal nuclease digestion (MNase-seq)



X-ChIP


High resolution mapping of protein-chromatin associated proteins



ORGANIC


Occupant regions of genomes from affinity-purified or affinity-enriched chromatin (ORGANIC)





Illumina sequencing workflow

- Library Construction
- Cluster Formation
- Sequencing
- Data Analysis



Fragmentation

- Mechanical shearing:

- BioRuptor
- Covaris

DNA, RNA


- Enzymatic:

- Fragmentase, RNAse3


DNA, RNA

- Chemical: Mg^{2+} , Zn^{2+}

→ RNA



DNA library construction



Fragmented DNA

↓ End Repair




Blunt End Fragments

↓ "A" Tailing




Single Overhang Fragments

↓ Adapter Ligation



DNA Fragments
with Adapter Ends

Enrichment of library fragments




Illumina Sequencing Technology


Sequencing By Synthesis (SBS) Technology

DNA
(0.1-1.0 ug)


**Library
preparation**




Cluster generation



Sequencing




TruSeq Chemistry: Flow Cell





GCTGAAAC
ATTCCCT
TCTGGGA
GAAATTAT
TGTTGA
AAGGAG
TTTGGG
CGCCAG
TCCCGC
AATTGC
TCTCC
AAGGCT
AATTGA
GCACAA
ATACCA
GCTTTT
TTTATZ


Sequencing





Sequencing



Patterned Flowcell




Hiseq 3000: 478 million nanowells per lane




CONCEPTUAL CHALLENGE— BEATING POISSON

Amplification Phase

Polyclonal (non-PF) Pads



Maximizing Well Occupancy and Monoclonality





Poisson statistics limit max
monoclonal occupancy < 40%

Polyclonality rises as occupancy
increases


SIMULTANEOUS SEEDING AND AMPLIFICATION

Maximizing Well Occupancy and Monoclonality



Amplification occurs at rate \gg faster
than seeding rate

Templates excluded from occupied
wells



What will go wrong ?

- cluster identification
- bubbles
- synthesis errors:

```
ClusterCluster  
ClustsrCluster  
ClusterCluster  
ClusterCluster  
Cl1sterCluster
```

What will go wrong ?

➤ synthesis errors:


ClusterCluster
ClustsrCluster
ClusterCluster
ClusterCluster
Cl1sterCluster

ClsterClusterC
ClusterCluster
ClusterCluster
Cl1usterCluste
ClusterCluster


Phasing & Pre-Phasing
problems

Illumina SAV viewer


TTGAGATT
TATGAAGC
TAAGTCTC
TAGAAGCC




base composition




fluorescence intensity




amplicon mix




amplicon




amplicon mix Q30




A Review of Sequencing Duplicate Types



**“If you can put adapters on it,
we can sequence it!”**



Know your sample



No need to be scared of HTS

UC Davis Center for Plant Diversity/Herbarium

- The Herbarium archives contain over 300,000 dried specimens.
- Search for **Grapevine Red Blotch-Associated Virus**
- Virus traces found by PCR



Maher Al Rwahnih
UCD Plant Foundation
Plant Services

Studying historic Bean varieties from herbarium samples

- GBS (Genotyping-By-Sequencing)
- 60 year old herbarium samples



Sarah Dohle,
Gepts Lab

Quantitation & QC methods

➤ Intercalating dye methods (PicoGreen, Qubit, etc.):

Specific to dsDNA, accurate at low levels of DNA

Great for pooling of indexed libraries to be sequenced in one lane

Requires standard curve generation, many accurate pipetting steps

➤ Bioanalyzer:

Quantitation is good for rough estimate

Invaluable for library QC

High-sensitivity DNA chip allows quantitation of low DNA levels

➤ qPCR

Most accurate quantitation method

More labor-intensive


Must be compared to a control

Optional: PCR-free libraries



- PCR-free library:
 - if concentration allows
 - Reduction of PCR bias against e.g. GC rich or AT rich regions, especially for metagenomic samples

OR

- Library enrichment by PCR:
 - Ideal combination: high input and low cycle number; low-bias polymerase




Library QC by Bioanalyzer




Predominant species of appropriate MW
Minimal primer dimer or adapter dimers
Minimal higher MW material


Library QC by Bioanalyzer



Beautiful




100% Adapters



Beautiful

Library QC




Examples for successful libraries



Adapter contamination
at ~125 bp

RNA-seq targeted sequencing:

- Capture-seq (Mercer et al. 2014)
- Nimblegen and Illumina
- Low quality DNA (FFPE)
- Lower read numbers 10 million reads
- Targeting lowly expressed genes.





<http://pacificbiosciences.com>


THIRD GENERATION DNA SEQUENCING



Single Molecule Real Time (SMRT™) sequencing
Sequencing of single DNA molecule by single
polymerase

Very long reads: average reads over 8 kb, up to 30 kb
High error rate (~13%).


Complementary to short accurate reads of Illumina



Third Generation Sequencing : Single Molecule Sequencing

Pacific Biosciences



70 nm aperture
“Zero Mode
Waveguide”




4 nucleotides with different
fluorescent dye simultaneous
present

2-3 nucleotides/sec
2-3 Kb (up to 50) read length
6 TB data in 30 minutes

laser damages polymerase






First Sequencing of CGG-repeat Alleles in Human Fragile X Syndrome using PacBio RS Sequencer

Paul Hagerman, Biochemistry and Molecular Medicine, SOM.

- Single-molecule sequencing of pure CGG array,
 - first for disease-relevant allele. Loomis *et al.* (2012) *Genome Research*.
 - applicable to many other tandem repeat disorders.
- Direct genomic DNA sequencing of methyl groups,
 - direct epigenetic sequencing (paper under review).
- Discovered 100% bias toward methylation of 20 CGG-repeat allele in female,
 - first methylated DNA sequencing in human



Future's so bright



TTGAGATT
TATGAAGC
TAAATCTC
TACCCCT
GCTGAAB
ATTCCCT
TCTGGGA
GAAATTAT
TGTTGA
AAGGAG
TTTGGG
CGCCAG
TCCCGC
AATTGC
TCTCC
AAGGCT
AATTGA
GCACAA
ATACCA
GCTTTT
TTTATAT



Thank you!

